

BRIDGING THE GAP

Impact and Correction of
Camera Noise for Computational
Microscopy including Precision
Localization Nanoscopy

Keith Bennett, Ph.D., Stephanie Fullerton, Ph.D.,
Eiji Toda, Hiroyuki Kawai, Teruo Takahashi

ADAPTED FROM PRESENTATION GIVEN AT METHODS AND APPLICATIONS OF FLUORESCENCE,
GENOA, ITALY
SEPTEMBER 10, 2013

HAMAMATSU PHOTONICS K.K.
System Division


Cameras are **NOT** perfect!

Why is a
camera manufacturer
proclaiming
that
cameras are not perfect?

Because NO camera is perfect
&
Because understanding why
matters to your science

WHAT IS THE GAP?

The difference between the performance of an actual camera and a theoretically perfect camera



THE AMOUNT OF THE GAP DEPENDS ON:

1. Sensor technology
 - CCD
 - EMCCD
 - sCMOS
2. Camera specifications
 - Quantum Efficiency
 - Camera Noise
 - Read noise
 - Excess Noise
 - Photo-response non-uniformity (PRNU)
3. Input photon level
 - Ultra low light
 - Low Light
 - Intermediate
 - High

UNDERSTANDING **WHY** THERE IS A GAP ENABLES:

- Appropriate camera selection
- Optimized camera usage
- Optimized experimental design
- **More reliable data analysis**




Better Results

THE REAL CAMERAS

CCD


- Well established technology
- All electron to digital conversion done in one chain
- Limited speed
- Moderate read noise
- Very low dark current
- High QE
- Best pixel response uniformity



THE REAL CAMERAS

EMCCD


- Back-thinned for increased QE
- High voltage gain register on sensor to achieve on-chip amplification
- All electron to digital conversion through one chain (either for EM or no EM)
- Read noise is low due to gain
- Stochastic EM amplification adds excess noise and long tail



THE REAL CAMERAS

CMOS

- Newest technology
- Every pixel and column has own amplifier
- Very low mean rms read noise
- **Pixel dependent read noise**
- Fastest speeds and largest field of view
- FPGA processing achieves excellent response uniformity (low PRNU)



SEEING THE PREDICTED GAP

Single Pixel Noise & SNR

Fixed Pattern Noise &
Image SNR

(from **specs**)

THE PERFECT CAMERA


100% QE { Every photon is converted into one electron

0 e- read noise { Every electron is digitized exactly as expected every time


0% fixed pattern noise { Every pixel and amplifier perform identically and predictably

In a perfect camera, the SNR of a single pixel is limited only by the physics of photon statistics... i.e. shot noise.

$$SNR = \sqrt{S}$$




EFFECT OF QE ON THE GAP



A reduction in QE reduces SNR at all light levels

RELATIVE SNR (rSNR) PLOTS CLEARLY SHOW THE GAP



rSNR is the SNR for a camera plotted relative to the perfect camera

rSNR shows differences among cameras over full range of signal level

{ All SNR graphs in this talk will be presented as rSNR }

THE SIMPLE SIGNAL TO NOISE RATIO (SNR)

$$SNR = \frac{QE * S}{\sqrt{F_n^2 * QE * (S + I_b) + (N_r / M)^2}}$$

QE: Quantum Efficiency

S: Input Signal Photon Number (photon/pixel)

F: Noise Factor (= 1 for CCD/sCMOS and √2 for EM-CCD)

N_r : Readout Noise

M: EM Gain (=1 for CCD / CMOS)

I_b : Background

“Changing the Game”

EMCCDs: EXCESS NOISE IS THE REASON FOR THE GAP

SNR for CCD / CMOS



$$SNR = \frac{QE \times P}{\sqrt{QE \times P}}$$

$$= \sqrt{QE \times P}$$

SNR for EM-CCD



$$SNR = \frac{M \times QE \times P}{F_n \times M \times \sqrt{QE \times P}} = \sqrt{\frac{QE \times P}{F_n^2}}$$

$$= \sqrt{QE_{eff} \times P}$$

$$QE_{eff} = \frac{QE}{F_n^2} = \frac{QE}{2}$$

QE: Quantum Efficiency,
 P: Input Signal Photon Number,
 M: EM Gain
 F_n: Noise Factor
 (assumes dark current and read noise are negligible)


THE REAL CAMERAS

	CCD	EMCCD	CMOS
Sensor Type	Charged Coupled Device Interline	Electron Multiplying CCD Back-thinned	Complimentary Metal Oxide Sensor
Camera Name	ORCA-R2	ImagEM x2	ORCA Flash4.0 V2
Pixel Number	1024 x 1344	512 x 512	2048 x 2048
Pixel Size	6.45 μm x 6.45 μm	16 μm x 16 μm	6.5 μm x 6.5 μm
QE (@650 nm)	58 %	90 %	72 %
Frame Rate	18 fps / 8 fps	70 fps	100 fps / 30 fps
Relative read noise (N_r/M), single-frame rms	10 e ⁻ / 6 e ⁻	< 0.2 e ⁻ (M = 200)	1.9 e ⁻ / 1.3 e ⁻
Noise Factor (F_n)	1	$\sqrt{2}$ @ M>10	1

MIND THE GAP: PREDICTED PIXEL rSNR PERFORMANCE FOR THE MOST COMMON CAMERAS


1. { A camera with the highest SNR at the lowest light level may not be the best at higher light levels

2. { The SNR of an EMCCD above 1 electron/pixel is comparable to a camera with $QE_{eff} = QE/2$ due to excess noise from EM gain.



BEYOND THE SWEET SPOT: THE GAP EXPANDS AT HIGH LIGHT IF PRNU IS NOT CORRECTED

Single Frame rSNR



- PRNU reduces SNR at high light
- Cannot be subtracted from image
- “Raw” PRNU varies by sensor
- Can be corrected in camera to varying degrees

Model:

- QE: 70%
- Noise Factor (F_n): 1
- Read Noise 3 photons rms
- PRNU (σ_s): 1%

Image
SNR

$$SNR = \frac{QE * S}{\sqrt{\underbrace{F_n^2 * QE * (S + I_b)}_{Shot\ Noise} + \underbrace{\left(\frac{Nr}{M}\right)^2}_{Read\ Noise} + \underbrace{[\sigma_s * QE * (S + I_b)]^2}_{Fixed\ Pattern\ Noise}}}$$

All SNR curves will be rSNR @ $\lambda=650$ nm



MEASURING THE REAL GAP

An in-depth look at noise in
CCD, EMCCD and CMOS
cameras

ORCA-R2 INTERLINE CCD: PREDICTABLE AND ROBUST


1. { PRNU is insignificant

2. { Read noise histogram has single Gaussian distribution



Read Noise
(N_r/QE)


Shot Noise
& QE




EMCCD: SOME SURPRISING RESULTS

- Thickness variations from back-thinning process causes **spectrally-dependent** PRNU

 - Cannot be removed during manufacturing
 - Must be calibrated by users for *their specific spectrum*.
 - Individual pixel map required for correction




- The Gap for EMCCD in CCD mode becomes very wide due to PRNU



COMPLEX BEHAVIOR: A CLOSER LOOK AT EMCCD SNR WITH HIGH AND LOW GAIN

Complex Behavior


- Excess noise (eQE)
- PRNU
- Saturation
- High read noise (34 e- @ M=5, 70 fps)
- Gain hard to measure




EMCCD GAIN CAUSES UNEVEN PROBABILITY DISTRIBUTIONS

In simulated probability distribution functions for EMCCD, the output at high gain is **not** Poisson due to the electron multiplication process!


2 Photon Average Input
Gain = 200




10 Photon Average Input
Gain = 200



ORCA-FLASH4.0 V2 (sCMOS): A VERY COMFORTABLE SWEET SPOT




NOT ALL CAMERAS ARE CREATED EQUAL: FLASH4.0 SWEET SPOT BROUGHT TO YOU BY HAMAMATSU CAMERA ENGINEERS




Signal amplified and digitized in column-parallel ADC.
FPGA provides offset and gain correction to the raw digitized signal.

sCMOS: PIXEL-DEPENDENT READ NOISE

Rms read noise matches single frame rSNR.



Single Frame Dark Histogram




Does not fit a Gaussian distribution, i.e. is not completely modeled by a single "read noise."

BRIDGING THE GAP



Using knowledge of camera noise to get highest camera performance provides improved Precision Localization Results

WHAT IS
MOST
IMPORTANT?



throughput

field of view

sample contrast

frame rate

RESOLUTION

ACCURACY

SAMPLE **BRIGHTNESS**

distance measurement

background

minimum bleaching rate

TWO PHASES OF PRECISION LOCALIZATION MICROSCOPY

1. Collect Image Data



- Prepare Sample
- Minimize Background
- Optimize Optical System
- Consider Camera Induced Noise

2. Reconstruct Superresolution Image



- Calibrate Camera
- Implement Noise Corrections
- Apply Statistical Algorithms

STANDARD PRACTICE IS NOT THE BEST PRACTICE: USING EMCCD WITH GAIN YIELDS LEAST ACCURATE RESULTS

Mean photon count	Ultimate accuracy limit (nm)	Conventional EMCCD accuracy limit (nm)	UAIM at 900× accuracy limit (nm)	UAIM at 4500× accuracy limit (nm)	CCD accuracy limit ^a (nm)
200	4.84	9.71 [100.6%]	5.41 [11.8%]	5.13 [6.0%]	7.94 [64.0%]
400	3.42	6.94 [102.9%]	3.92 [14.6%]	3.63 [6.1%]	5.26 [53.8%]
800	2.42	4.93 [103.7%]	2.85 [17.8%]	2.57 [6.2%]	3.45 [42.6%]
1600	1.71	3.49 [104.1%]	2.08 [21.6%]	1.82 [6.4%]	2.32 [35.7%]
3200	1.21	2.48 [105.0%]	1.52 [25.6%]	1.30 [7.4%]	1.57 [29.8%]

^aComputed at near-optimal magnification (i.e., magnification that yields approximately the best localization accuracy limit) of 128.6×, 185.7×, 185.7×, 185.7×, and 242.9× for mean photon count of 200, 400, 800, 1600, and 3200, respectively.

CCD QE: 100%, read noise = 1.8 ph, **no background; No fixed pattern noise.**


... “the fact that the noise coefficient approaches 1 with increasing photon count demonstrates the suitability of the CCD (CMOS) detector when enough light is available....”

... “even with a low readout noise of $\sigma=2$ electrons, the CCD detector is unsuitable for extreme low light imaging and implementing UAIM.”

Adapted from: J. Chao et al (Ober Lab), Nat. Meth10, 2013) doi:10.1038/nmeth.2396

<http://www.wardoberlab.com/>

COMPENSATING READ NOISE VARIATION



Incorporating pixel-specific read noise into the Maximum Likelihood Probability Model eliminates and narrows the asymmetric distribution of localized molecules caused by higher read noise pixels.

Courtesy F. Huang, Bewersdorf Lab

Changing the Game

SELECTING AND USING CAMERAS: CASE STUDIES

CCD

Subnanometre single-molecule localization, registration and distance measurements

Alexandros Pertsinidis^{1,2}, Yunxiang Zhang^{1,2} & Steven Chu^{1,2,3,4,†}

EMCCD

Ultrahigh accuracy imaging modality for super-localization microscopy

Jerry Chao¹⁻³, Sripad Ram¹⁻³, E Sally Ward² & Raimund J Ober^{1,2}

sCMOS

Video-rate nanoscopy using sCMOS camera-specific single-molecule localization algorithms

Fang Huang¹, Tobias M P Hartwich^{1-3,9}, Felix E Rivera-Molina^{1,9}, Yu Lin^{4,5}, Whitney C Duim¹, Jane J Long⁶, Pradeep D Uchil⁷, Jordan R Myers¹, Michelle A Baird⁸, Walther Mothes⁷, Michael W Davidson⁸, Derek Toomre¹ & Joerg Bewersdorf^{1,4,5}

Subnanometre single-molecule localization, registration and distance measurements


Alexandros Pertsinidis^{1,2}, Yunxiang Zhang^{1,2} & Steven Chu^{1,2,3,4,†}

Results	Accurate measurement of the <i>distance</i> between two fluorophores of different colors. $\sigma_{\text{distance}} \sim 0.77$ nm using a dichroic beamsplitter to direct each color of light to separate halves of the CCD camera.	
Camera Correction	Measured PRNU maps for each color. Improved localization relative accuracy by ~ 2 –4 nm.	
Details	Speed: 5 – 50 s / measurement Light: $\sim 4,000$ – 10,000 ph/ mol/frame $\sim 10^5$ ph / mol before bleaching	Imaging: Simultaneous 2 color Camera: Back-thinned EM-CCD, gain off

Nature (2010) | doi:10.1038/nature09163

Ultrahigh accuracy imaging modality for super-localization microscopy

Jerry Chao¹⁻³, Sripad Ram¹⁻³, E Sally Ward² &
Raimund J Ober^{1,2}



Cholera toxin B subunit

scale bar: 1 μm

Results

Localization Microscopy with Minimal Bleaching. Plasma membrane dynamics for > 60 s (594 frames). 40% better localization precision than “conventional” EMCCD localization

Camera Correction

Implemented detailed statistical EM noise model into maximum likelihood reconstruction probability model.

Details

Speed: ~60s / reconstructed image **Mag:** 630X
Light: ~100 photons /molecule frame **Camera:** EM-CCD, Gain ~1000

Courtesy of J. Chao et al (Ober Lab)


Adapted from Nat Meth (2013) doi:10.1038/nmeth.2396

Video-rate nanoscopy using sCMOS camera-specific single-molecule localization algorithms

Fang Huang¹, Tobias M P Hartwich^{1-3,9}, Felix E Rivera-Molina^{1,9}, Yu Lin^{4,5}, Whitney C Duim¹, Jane J Long⁶, Pradeep D Uchil⁷, Jordan R Myers¹, Michelle A Baird⁸, Walther Mothes⁷, Michael W Davidson⁸, Derek Toomre¹ & Joerg Bewersdorf^{1,4,5}

Results

32 reconstructed frames / sec ($6.6 \times 6.6 \mu\text{m}^2$) field of view; fixed and living cells showed cellular dynamics not visible in reconstructions using longer data collection times.



Camera Correction

Implemented pixel-specific read noise into probability model for MLEM.


Details

Speed: 0.03 s/ reconstructed image
Light: ~3,000 photons /mol/ frame

Mag: 60X,
Camera: sCMOS, 3200 fps

Courtesy of F. Huang, (Bewersdorf Lab)

Localization Precision “conventional” EMCCD vs. sCMOS




Courtesy of F. Huang. Bewersdorf Lab, Yale

Adapted from F. Huang *et al.*, Nature Methods 10(7): 653-658 (2013)

MINIMIZING THE GAP: MATCHING THE CAMERA TO YOUR NEEDS

Higher Accuracy
Better Resolution
Lower Sample Contrast



HAVE YOU DONE A **GAP** ANALYSIS?

1. How much light do I have?

The relative performance of CCD, back-thinned CCD, EMCCD and sCMOS cameras is light-level dependent.

2. Do I know my camera's strengths and weaknesses?

No camera is perfect; proper use is required for the best results and to avoid errors

3. What is the goal of my experiment?

The most appropriate choice of camera depends upon the specific super resolution / localization experiment

Acknowledgements

Prof. Zhen-li Huang, Huanzhong University of Science and Technology

F. Long et al, OPTICS EXPRESS 17741 (2012)

Prof. Joerg Bewersdorf, Yale University

F. Huang *et al.*, Nature Methods 10(7): 653-658 (2013) See a movie of 32 fps dynamics in the supplementary materials of the article by Huang et al.

Prof. Raimund Ober, Texas Southwestern University

J. Chao et al, Nat Meth (2013) doi:10.1038/nmeth.2396

Hamamatsu

Hiroyuki Kawai: camera measurements

Teruo Takahashi: simulations

Stephanie Fullerton: presentation preparation

Keith Bennett

kbennett@hamamatsu.com