

Photon number resolving capability of qCMOS camera for Raman spectroscopy and imaging

Kenichiro Ikemura, Kunihiro Tsuchiya, Taiki Miura
Systems Division, Hamamatsu Photonics K.K., Japan

HAMAMATSU
PHOTON IS OUR BUSINESS

Introduction of qCMOS

Introducing the new qCMOS

ORCA-Quest




- Unprecedented Performance
- Ultimate sensitivity
- Photon number resolving capability

LOW READOUT NOISE 0.27 ELECTRONS RMS <small>Ultra-Quiet Scan</small>	HIGH QE 90% @475 nm <small>Back-illuminated qCMOS</small>	HIGH RESOLUTION 4096 × 2304 <small>9.4 Megapixels</small>	HIGH SPEED 120 FRAMES/S <small>@4096 × 2304 (16 BIT) STANDARD SCAN</small>
-----------------------------------------------------------------------------------	------------------------------------------------------------------------	------------------------------------------------------------------------	-----------------------------------------------------------------------------------------

What does "Photon number resolving" mean?


- Photon number resolving is unique and quite different from single photon counting.
- The existing photon counting cameras like EM-CCD can never do photon number resolving because of their excess multiplication noise.
- The following graph shows the probability distribution of the observed photoelectrons with a mean of $N = 3$ photoelectrons, and three different values of the readout noise. The smaller the readout noise is the deeper are the valleys between the photoelectron peaks. If the CMOS sensor can reduce the readout noise, it enables photon number resolving.



qCMOS capability for Raman imaging

Raman imaging method qCMOS can be beneficial

	1. Point scan method	2. Line scan method	3. Area scan method
Detection area	Small	Middle	Large
Scan speed	Slow	Middle	Fast
Spectral resolution	High	High	Low
Spatial resolution	Depends on scan pitch	Depend on scan pitch	Depend on camera
Focus	Good	Depend on sample	Depend on sample
Detection	Spectrometer + line sensor	Spectrometer + camera	Band-pass filter + camera
Products	Many	Some	Very few



1. Point scan method : **NG**
qCMOS has too many vertical pixels with small pixel size compared to common sensors for Raman. Digital binning of CMOS camera accumulates readout noise of each pixel, which might lead to lower sensitivity, low throughput.


2. Line scan method : **OK**
The method does not use binning function, so we can enjoy benefit from qCMOS such as pixel size, pixel number, low noise, high QE etc.

3. Area scan method : **OK**
Same as the line scan method. This method could be evaluated from ordinal SN simulation.

We want to evaluate the benefit quantitatively by using Quest on "line scan method". Especially, "spectral SN" is practical for the evaluation.

Assumed line scan type Raman imaging systems


- We tested what the effect can be brought by replacing EMCCD, a common camera for Raman imaging, to qCMOS for line scan type Raman imaging system.
- 1st assuming case is **retrofit system** where there is no change of the optics and spectrometer and only a camera is replaced.
- 2nd assuming case is **customized system** where the optics and spectrometer is optimized for qCMOS.




Evaluation of SN, SN_{sp}

How to evaluate SN (Signal-to-noise ratio)

- Create a data set of 30 images
- Create a spectrum from a single line in the center of the image
- We evaluate SN and spectral SN as shown below. (I: Intensity, σ : standard deviation)
 - $SN = \frac{I_{peak}}{\sigma_{peak}}$ (Conventional SN)
 - $SN_{sp} = \frac{I_{peak}}{\sigma_{Base}}$ (Spectrum SN)




Spectrum comparison under condition of same PN (photon number) / pixel




*1) The term "PNR" means Photon number resolving. PNR mode outputs the number of photoelectrons from Digital output.
*2) The term "4Bin" means treating 16 adjacent pixels as one big pixel (arranged in 4x4 pixels)

SN comparison under same PN/pix condition




SN_{sp} comparison under same PN/pix condition



Spectral dynamic range


An important point of spectral analysis is how many Raman bands in wide intensity range can be detected.

- Minimum detection limit is determined by spectrum SN, SN_{sp} . Assume that the Raman signal can be detected when $SN_{sp} > 2$.
- Maximum detection limit can be adjusted by high intensity Raman band and maximum AD counts of a detector




Customized system

EM-CCD vs qCMOS vs qCMOS-PNR



Retrofit system


EM-CCD vs qCMOS-PNR (4Bin)



- The spectral dynamic range of qCMOS is larger than that of EM-CCD.
- The spectral dynamic range of EM-CCD is deteriorated by EM gain amplification.*
- The spectral dynamic range is shown in the table below. *In this experiment, EM gain is 640x.

Dynamic range	EM-CCD	qCMOS	qCMOS-PNR	qCMOS-PNR (4Bin)
PN	2.3~570	1.1~10K	0.8~10K	1.2~160K
AD count	240~60K	7~60K	5~60K	7~970K



- The figure attached below shows the spectral dynamic range when measurement conditions are set to the maximum AD count for a sample with no background.



Quest is higher dynamic range than EM-CCD >> As a model of customized system
Quest-PNR (4Bin) is higher dynamic range than EM-CCD because of PNR-binning effect. >> As a model of retrofit system

Samples with background emission

- For the SN evaluation, we used acetone which do not have background emission such as fluorescence.
- In this section, we measured Raman spectrum of Lycopene contained in tomato juice, which does have fluorescence around Raman resonance wavelength.



- The Raman spectrum by qCMOS cameras is better than EM-CCD.
- EM-CCD amplifies both of signal and background. Therefore, EM-CCD noise (especially by a contribution of excess noise) would be higher than qCMOS camera noise.

Conclusions

Table: Consideration of line scan type Raman imaging system performance of replacing EM-CCD to qCMOS

Contents	Retrofit	Retrofit (w binning&PNR)	Custom design
Mechanics			
Lens	—	—	Use a proper designed lens
Spectrometer	—	—	Use a short focal length or low ruling number of grating
Input split	—	—	Depends on measurement design
Adaptor	Custom made	Custom made	Custom made
Spatial resolution	Better	Similar	Similar
Detectable light intensity			
SN	Worse	Better or Similar	Better
Speed	Worse	Better or Similar	Better
Spectral resolution	Similar	Similar	Similar
Spectral dynamic Range	Better	Better	Better
FOV range	Better or Similar	Better or Similar	Better or Similar

- For retrofit system, the performance of qCMOS-PNR (4Bin) is better than EM-CCD.
- For customized system, the performance of qCMOS and qCMOS-PNR is better than EM-CCD.